



i-gel[®] O₂ Resus Pack

Three things in one pack, one less thing to think about



Airway Management ▪ Airway Devices



In emergency medicine, you need equipment that's easy, rapid and reliable to use. The i-gel® O₂ Resus Pack contains everything you need to prepare, insert and secure the i-gel® O₂ quickly and efficiently: an i-gel® O₂ supraglottic airway, a sachet of lubricant and an airway support strap.

The i-gel® O₂ incorporates a supplementary oxygen port and an integral hook ring for securing of the airway support strap.

The i-gel® O₂ has been designed to facilitate ventilation as part of standard resuscitation protocols, such as those designated by the European Resuscitation Council (ERC) and the American Heart Association (AHA). However, the i-gel® O₂ incorporates a supplementary oxygen port, so that it can also be used for the delivery of passive oxygenation, or Passive Airway Management (PAM™), as part of an appropriate CardioCerebral Resuscitation (CCR) protocol.



Airway innovation since 2007

i-gel® is the innovative second generation supraglottic airway device from Intersurgical launched in 2007. The first major development since the laryngeal mask airway, the i-gel® changed the face of airway management and is now widely used in anesthesia and resuscitation around the world.



The non-inflating cuff

The i-gel® gets its name from the innovative soft, gel-like material from which it is made. It is the innovative application of this material that has enabled the development of a unique non-inflating cuff. This means there is

no need for cuff deflation prior to insertion and no cuff inflation after placement in order to secure a seal, which shortens and simplifies the preparation and insertion procedure.

Insertion

The i-gel® O₂ is incredibly easy to use. Insertion is rapid and can normally be achieved in less than 5 seconds, by a proficient user. The simplicity of i-gel® O₂ means the insertion technique is instinctive, reliable and does not require or benefit from the insertion of fingers or thumbs into the patient's mouth.

Fixation

The i-gel® O₂ Resus Pack includes a specially designed airway support strap for securing the i-gel® O₂ in position. Fixation is easily achieved by sliding the strap under the patient's neck and securing it in place with the hook ring on the collar of the i-gel® O₂. This makes it ideal for use where adhesive tape is unsuitable, such as on patients with facial hair.



Visit the i-gel® website
www.i-gel.com



i-gel[®] O₂ Resus Pack

15mm connector

Reliable connection to any standard catheter mount or connection

Proximal end of gastric channel

Color coded hook ring

To secure the i-gel[®] O₂ in position with the airway support strap

Clearly displayed product information

For quick easy reference. Includes confirmation of size and weight guidance

Position guide

(adult sizes only)

Supplementary oxygen port

For the administration of passive oxygenation as a component of cardiocerebral resuscitation (CCR)

Gastric channel

Enhances patient safety by providing a mechanism for the management of regurgitant fluid

Integral bite block

Reduces the possibility of airway channel occlusion

Buccal cavity stabilizer

Aids insertion and eliminates the potential for rotation

Epiglottic rest

Reduces the possibility of epiglottis 'down folding' and airway obstruction

Non-inflating cuff

Eliminates the need for cuff inflation after insertion, allowing easy and rapid insertion

Distal end of gastric channel



The i-gel[®] O₂ Resus Pack is available in 3 adult sizes and includes:




- i-gel[®] O₂ supraglottic airway
- Sachet of lubricant – for quick and easy lubrication of the i-gel[®] O₂ prior to insertion
- Airway support strap – to secure the i-gel[®] O₂ in position



Ordering information

The i-gel® O₂ Resus Pack contains an i-gel® O₂, a sachet of lubricant and an airway support strap all packaged in a specially designed clear, sterile, rigid pack. This unique design ensures protection of the contents, both in transit and in storage. The i-gel® O₂ has a color coded hook ring to allow quick and easy identification of size. The i-gel® O₂ Resus Pack is available in three adult sizes.



Code	Description	Weight	Qty.
8705030	i-gel® O ₂ Resus Pack, large adult – includes a size 5 i-gel® O ₂ with orange hook ring, sachet of lubricant and airway support strap	90+kg (200+lbs)	6 
8704030	i-gel® O ₂ Resus Pack, medium adult – includes a size 4 i-gel® O ₂ with green hook ring, sachet of lubricant and airway support strap	50–90kg (110–200lbs)	6 
8703030	i-gel® O ₂ Resus Pack, small adult – includes a size 3 i-gel® O ₂ with yellow hook ring, sachet of lubricant and airway support strap	30–60kg (65–130lbs)	6 

[Make an inquiry](#)
[Find out more](#)
[Watch the video](#)

 Sterile

US Lit No. 100.4 • 03.23