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en INSTRUCTIONS FOR USE

Intended use

The PEEP (Positive End Expiratory Pressure) valve is intended for use in respiratory systems to maintain positive end expiratory pressure in the airways.

Restrictions of use

- The device must be used by qualified and trained medical/nursing staff aware of the possible side effects of positive pressure ventilation.
- To be used for administration of air and oxygen.
- If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.

Contraindications

Non-invasive positive pressure ventilation is contraindicated in the following cases:

- Coma
- uncooperative patient
- cardiac arrest
- hemodynamic instability
- recent oesophageal and gastro-surgical operations
- Heavy bleeding of the upper digestive tract
- obstruction of the upper airways
- Pneumothorax

Complications/side effects

The efficacy of positive pressure ventilation is considerably affected by the PEEP values set. Therefore, the medical staff is advised to carefully evaluate the pressure level most suited to the clinical condition of the patient. Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar overdistension.

Instructions for use

Remove the device from the package. Check correct movement of the knurled cap of the adjustable PEEP valve. Check the movement of the internal plate before connecting the device. Connect the PEEP valve via ISO standard connections on relevant equipment such as hood or mask; the device can be used with CPAP breathing systems and CPAP hoods. Adjustable PEEP: set the desired pressure by acting on the cap (clockwise to increase the PEEP, anticlockwise to decrease it). Possible adjustments from 0 - 20 cm H2O (mbar). Pre-calibrated PEEP: Always check that the pressure indicated on the device is as desired.

The device is now ready for operation. Proceed with delivery of the ventilation gases.

Warnings and precautions

⊗ Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

- In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
- The measurement (approximate) that the device provides is not significantly affected by variation of the ventilation flows.
- It is recommended to check the PEEP with a pressure gauge.
- Monitor the clinical parameters of the patient during the treatment.
- Use only with devices that have connections in compliance with the ISO standards.
- Caution: The device contains metal; do not use in departments where the presence of metal is a source of risk for the safety of the patient or third parties.
- The device can be used for maximum 7 days. Depending on the clinical picture of the patient the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
- Non-sterile. Do not sterilize.
- Do not use on another patient.
- Expiry: 5 years provided that the packa- ging is undamaged and if stored in normal storage conditions [-20°C/+50°C].
- Dispose of the materials immediately after use CONFORMING TO the current laws and regulations.
- Use only with devices that have connections in compliance with the ISO 5356-1 standards.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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