



VENTUPLUS

- (en) ACCESSORY FOR HELMETS AND MASKS - VENTUPLUS
- (fr) ACCESSOIRE POUR HEAUMES ET MASQUES DE VENTILATION - VENTUPLUS
- (de) ZUBEHÖR FÜR VENTUPLUS BEATMUNGSELME UND - MASKEN
- (es) ACCESORIO PARA CASCOS RESPIRATORIOS Y MASCARILLAS - VENTUPLUS
- (pt) ACESSÓRIO PARA HELMET (CAPACETE) E MÁSCARAS - VENTUPLUS
- (it) ACCESSORIO PER CASCHI E MASCHERE - VENTUPLUS
- (nl) ACCESSOIRE VOOR HELMEN EN MASKERS - VENTUPLUS
- (no) TILBEHØR FOR HJELMER OG MASKER - VENTUPLUS
- (fi) LISÄLAITE KYPÄRIIN JA MASKEIHIN - VENTUPLUS
- (sv) TILLBEHÖR FÖR HJÄLMAR OCH MASKER - VENTUPLUS
- (da) TILBEHØR TIL HJELME OG MASKER - VENTUPLUS
- (el) ΒΟΗΘΗΤΙΚΟ ΕΞΑΡΤΗΜΑ ΓΙΑ ΚΑΣΚΕΣ ΚΑΙ ΜΑΣΚΕΣ ΑΕΡΙΣΜΟΥ - VENTUPLUS
- (pl) AKCESORIA VENTUPLUS DO HELMÓW I MASEK
- (cs) PŘÍSLUŠENSTVÍ PRO KUKLY A MASKY - VENTUPLUS
- (hu) KIEGÉSZÍTŐ A LÉLEGEZTETŐ SISAKOKHOZ ÉS MASZKOKHOZ - VENTUPLUS
- (sl) DODATEK ZA ČELADE IN MASKE - VENTUPLUS
- (et) KIIVRITE JA MASKIDE TARIK - VENTUPLUS
- (ro) ACCESORIU PENTRU CĂȘTI ȘI MĂȘTI - VENTUPLUS
- (sk) PŘÍSLUŠENSTVO PRE PRILBY A MASKY - VENTUPLUS
- (hr) DODATAK ZA KACIGE I MASKE - VENTUPLUS
- (tr) KASKLAR VE MASKELER İÇİN AKSESUAR - VENTUPLUS
- (ar) ملحقات الخوذات والأقنعة - فينتوبلاس
- (he) VENTUPLUS אביזר לקסדות ומסכות
- (ua) ПРИЛАДДЯ ДЛЯ ШОЛОМІВ ТА МАСОК - VENTUPLUS
- (sr) DODATAK ZA KACIGE I MASKE - VENTUPLUS

Instructions for use / Mode d'emploi / Gebrauchsanleitung / Instrucciones de uso / Instruções de utilização / Istruzioni per l'uso / Gebruiksaanwijzing / Brukerveiledning / Käyttöohjeet / Bruksanvisning / Brugsanvisning / Οδηγίες χρήσης / Instrukcja stosowania / Návod k použití / Használati útmutató / Navodila za uporabo / Kasutusjuhend / Instrucțiuni de utilizare / Návod na použitie / Upute za upotrebu / Kullanna talimatlari / تعليمات الاستخدام / הוראות שימוש / Інструкції з використання / Uputstvo za upotrebu

NOTE: suggested configuration / REMARQUE: conseil de configuration / HINWEIS: Konfigurationsempfehlung / NOTA: Configuración sugerida /
 NOTA: sugestão de configuração / NOTA: suggerimento di configurazione / OPMERKING: aanbevolen configuratie / OBS: foreslått konfigurasjon
 / HUOMAUTUS: ehdotettu kokoonpano / ANMÄRKNING: Förslag på konfiguration / BEMÆRK: foreslået konfiguration / ΣΗΜΕΙΩΣΗ: συνιστώμενος
 τύπος / UWAGA: konfiguracja sugerowana / POZNÁMKA: doporučená konfigurace / MEGJEGYZÉS: javasolt beállítás / OPOMBA: predlagana
 konfiguracija / MÁRKUS: soovitatud seadistus / NOTÁ: recomandare de configurare / POZNÁMKA: odporúčané konfigurácie /
 NAPOMENA: predložena konfiguracija / NOT: Konfigurasyon önerisi / ملاحظة: اقتراح الإعداد / הערה - תצורה מומלצת /
 ПРИМІТКА. Рекомендована конфігурація / NAPOMENA: predložena konfiguracija

(en) CONFIGURATIONS:

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INTRODUCTION

The VENTUPLUS accessory is based on the Venturi tube operating principle: a constant oxygen flow blows through the tube which has a constriction (nozzle). The speed and pressure of the oxygen flow undergo a variation which allows aspirating environmental air from the area surrounding the Venturi tube. The total flow is a mixture of air and oxygen whose O₂% can be modified by using the additional oxygen source if present. The performance table attached to the product indicates how to obtain the combinations of O₂% and total patient flow at each PEEP level.

DESTINATION OF USE

The accessory is used in combination with helmets and masks for CPAP (Continuous Positive Airway Pressure) treatment. To be used in hospital environment when high flow of air/oxygen mixture at variable O₂% is required.

INDICATIONS FOR USE

To be used with helmets and masks for CPAP (Continuous Positive Airway Pressure) treatment. It can deliver a high flow of air/oxygen mixture according to the setting by healthcare personnel.

RESTRICTIONS OF USE

- The device must be used by qualified and trained medical/nursing staff.
- To be used for administration of air and oxygen.
- If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.
- To be used only with oxygen flowmeters compliant to the relevant standards and capable of delivering the flows reported on the product's performance table.
- Use only with oxygen flowmeters featuring 7mm diameter hose tail.
- To be connected with respiratory systems/patient interfaces featuring 22mm connections (compliant with ISO standards).

WHERE IT CAN BE USED

The VENTUPLUS accessory allows the delivery of high flow CPAP therapy using compressed oxygen sources (wall-mounted distribution systems or cylinders) available in hospital departments and ambulances.

HOW TO USE VENTUPLUS

1. Remove the VENTUPLUS accessory from the package. Connect the device to the breathing system/patient interface gas inlet connector. Verify safe connections.
2. Attach the oxygen input connector/s A and B to the flowmeter/s (see figure 1). Check proper connection before activating the gas supply.
3. Activate the flowmeter/s according to the performance table attached to the device.
4. To supply O₂% equal to 100% the Venturi valve has to be inactivated by sliding the transparent cover (see figure 2a) until the air inlet is occluded (lock to closed position

– see figure 2b).

5. Only when the air inlet/transparent cover is closed, the total flow to the patient is equal to the sum of the oxygen flows from the two flowmeters. To restore the operation of the Venturi valve, slide and lock the transparent cover to the opened position (figure 2a) and adjust therapy conditions according to the performance table.

NOTE: The cover must be fully engaged to the opened/closed positions.


Fig. 1


Fig. 2a

Fig. 2b

PERFORMANCE TABLE

Refer to the performance table attached to the device for information about the oxygen consumption and delivered total patient flow.

WARNINGS AND CAUTIONS

⚠ Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, the materials of the device could be damaged by disinfection methods, the product performances may be impaired developing a safety risk for patients.

1. Depending on the clinical condition of the patient the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
2. It is recommended to use the settings provided by Intersurgical within the performance table attached to the device. These settings ensure sufficient flows for proper CPAP therapy. Higher flow settings would generate unnecessary oxygen consumption and excessive noise.
3. Always check that the oxygen supply lines are safely

connected. In the event of oxygen tubes disconnection from the device, patient ventilation will fail and oxygen could be released into the environment.

4. In the events of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
5. It is recommended to monitor the clinical parameters of the patient during the therapy.
6. If the patient's conditions deteriorate or there is no improvement within the expected time, it is recommended to evaluate alternative ventilation techniques.
7. Do not use oxygen bubble humidifiers with the gas supply flowmeters.
8. Non-sterile. Do not sterilize.
9. Expiry: 5 years provided that the packaging is undamaged and the product is stored in normal storage conditions [-20°C / +50°C].
10. Dispose of the device immediately after use in compliance with the current laws and regulations.
11. Use only with devices that features connections in compliance with the ISO 5356-1 standards.

CONFIGURATIONS:

VENTUPLUS CP924 WITH DUAL OXYGEN FLOWMETER, SUGGESTED CONFIGURATION FOR O₂% ABOVE 40% (UP TO 100%) SEE FIG. 3

VENTUPLUS CP924 WITH TWO SINGLE OXYGEN FLOWMETERS, SUGGESTED CONFIGURATION FOR O₂% ABOVE 40% (UP TO 100%) SEE FIG. 4

VENTUPLUS CP924 WITH SINGLE OXYGEN FLOWMETER, SUGGESTED CONFIGURATION FOR O₂% EQUAL TO 40% SEE FIG. 5

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy