



Solus™ laryngeal mask airways

A comprehensive range of single use laryngeal mask airways



Airway Management ▪ Airway Devices



Solus™ laryngeal mask airways

A comprehensive range of single use laryngeal mask airways suitable for anesthesia and emergency medicine

Designed to meet the needs of a variety of clinical scenarios, all Solus laryngeal mask airways are single use, supplied sterile and do not contain natural rubber latex.

Low friction material
Atraumatic insertion & removal

Classic cuff shape
Renowned for reliability

Essential information
Clearly visible for quick reference

Integral inflation line
Reduces the risk of damage while in use








High quality valve
To ensure continual cuff integrity

Solus Standard

Using specially chosen materials, the Solus Standard is manufactured with the quality, flexibility and strength required to work effectively in a range of demanding clinical applications, including pediatric anesthesia. Available in a full size range, from adult to neonate.



Make an inquiry

Code	Description	Size	Weight	Qty.
8005000	Solus Standard, laryngeal mask airway	5 Large adult	70+kg	20 
8004000	Solus Standard, laryngeal mask airway	4 Medium adult	50–70kg	20 
8003000	Solus Standard, laryngeal mask airway	3 Small adult	30–50kg	20 
8025000	Solus Standard, laryngeal mask airway	2.5 Large pediatric	20–30kg	20 
8002000	Solus Standard, laryngeal mask airway	2 Small pediatric	10–20kg	20 
8015000	Solus Standard, laryngeal mask airway	1.5 Infant	15–10kg	20 
8001000	Solus Standard, laryngeal mask airway	1 Neonate	11<5kg	20 



Solus Satin

Solus Satin is our range of single use laryngeal mask airways featuring a softer airway tube providing increased flexibility.



Make an inquiry

Code	Description	Size		Weight	Qty.	
8005003	Solus Satin, laryngeal mask airway	5	Large adult	70+kg	20	
8004003	Solus Satin, laryngeal mask airway	4	Medium adult	50–70kg	20	
8003003	Solus Satin, laryngeal mask airway	3	Small adult	30–50kg	20	

Solus MRI-safe

A range of single use, sterile, laryngeal mask airways fitted with specially designed non-ferrous valves, guaranteed not to interfere with the magnet in an MRI scanner. The plastic valve has been selected to ensure full reliability throughout the shelf life of each device.



Make an inquiry

Code	Description	Size		Weight	Qty.	
8005002	Solus MRI-safe, laryngeal mask airway	5	Large adult	70+kg	20	
8004002	Solus MRI-safe, laryngeal mask airway	4	Medium adult	50–70kg	20	
8003002	Solus MRI-safe, laryngeal mask airway	3	Small adult	30–50kg	20	
8025002	Solus MRI-safe, laryngeal mask airway	2.5	Large pediatric	20–30kg	20	
8002002	Solus MRI-safe, laryngeal mask airway	2	Small pediatric	10–20kg	20	
8015002	Solus MRI-safe, laryngeal mask airway	1.5	Infant	15–10kg	20	
8001002	Solus MRI-safe, laryngeal mask airway	1	Neonate	11<5kg	20	



Solus Curve


Solus Curve is our laryngeal mask for those who prefer the insertion characteristics of a device with a curved airway tube.

[Make an inquiry](#)

Code	Description	Size		Weight	Qty.	
8005050	Solus Curve, laryngeal mask airway	5	Large adult	70+kg	20	
8004050	Solus Curve, laryngeal mask airway	4	Medium adult	50–70kg	20	
8003050	Solus Curve, laryngeal mask airway	3	Small adult	30–50kg	20	

Solus Flexible

Solus Flexible provides an ideal solution for shared airway procedures such as ENT, dental, oro-maxillary and eye surgery where the rigidity of airway devices can obscure surgical access. The wire-reinforced tube permits flexion without kinking and can be moved at any time mid-procedure without concern for the cessation of gas flow.

[Make an inquiry](#)

Code	Description	Size		Weight	Qty.	
8005001	Solus Flexible, wire-reinforced laryngeal mask airway	5	Large adult	70+kg	20	
8004001	Solus Flexible, wire-reinforced laryngeal mask airway	4	Medium adult	50–70kg	20	
8003001	Solus Flexible, wire-reinforced laryngeal mask airway	3	Small adult	30–50kg	20	
8025001	Solus Flexible, wire-reinforced laryngeal mask airway	2.5	Large pediatric	20–30kg	20	
8002001	Solus Flexible, wire-reinforced laryngeal mask airway	2	Small pediatric	10–20kg	20	

Sterile

US Lit No. 24.6 • 03.23