



i-view™ video laryngoscope

Video laryngoscopy wherever and whenever you intubate



Airway Management ▪ Airway Accessories



Introducing the i-view™ video laryngoscope

One size ▪ Single use ▪ Fully disposable

i-view is the new, single use, fully disposable* video laryngoscope from Intersurgical, providing the option of video laryngoscopy wherever you might need to intubate.

By incorporating a Macintosh blade, i-view can also be used for direct laryngoscopy. The technique for insertion is more familiar and instinctive than for devices with a hyper-angled blade.

i-view is ready to use seconds after removing from the packaging. The ergonomic design ensures i-view is easy to use, and the integral LCD screen provides an optimal view in a variety of light conditions.



i-view provides the opportunity of video laryngoscopy 'whenever and wherever you might need to intubate'

The availability of a video laryngoscope may be limited due to the cost implications of purchasing reusable devices for multiple sites. By combining all the advantages of a fully integrated video laryngoscope in a single use, disposable product, i-view provides a cost effective solution.

i-view is ideal for use in:

- pre-hospital
- emergency medicine
- resuscitation
- anesthesia
- difficult airway
- intensive care



Intensive care



Pre-hospital/resuscitation



Difficult airway

*The Intersurgical i-view is a single use, fully disposable device. Disposal should be in accordance with normal procedures for clinical waste, including landfill or incineration. Local protocols or regulations should be followed when disposing of i-view.



Features & benefits

Video laryngoscopy wherever and whenever you intubate

One size

allows you to hold less stock

Single use

reduces the risk of cross infection and avoids reprocessing costs

Fully disposable

no need to remove batteries prior to disposal*

Integral LCD screen

provides a view in a variety of light conditions

Ergonomic design

ensures i-view is comfortable and instinctive to use



On/Off switch

makes i-view ready to use seconds after removing from the packaging

Camera and LED light

creates optimal visualization



Macintosh blade

allows i-view to be used for direct as well as video laryngoscopy



*The Intersurgical i-view is a single use, fully disposable device. Disposal should be in accordance with normal procedures for clinical waste, including landfill or incineration. Local protocols or regulations should be followed when disposing of i-view.



i-view is individually packed in a specially designed rigid shell, sealed with a plastic film to ensure protection of the device, both in transit and in storage. Each pack can be opened by identifying the corner tab labeled with a "thumb symbol". By pulling back the top plastic film until the product can be easily removed.

The LOT number and expiry date of the device is located on the top edge above the LCD screen.



Code	Description	Qty.
8008000	i-view video laryngoscope	8

[Make an inquiry](#)[Find out more](#)[Watch the video](#)