




## Air-Guard™ Clear

Filter for use in respiratory systems and oxygen concentrators



**Airway Management** ▪ Breathing Filters, HMEs and HMEFs



## Air-Guard™ Clear

### For use in respiratory systems and oxygen concentrators


The Intersurgical Air-Guard Clear filter is a hydrophobic pleated membrane filter designed for the protection of oxygen concentrator machines and other respiratory equipment.

The Air-Guard Clear filter will function as a microbial filter when used to protect oxygen concentrator machines or when used in breathing systems to protect equipment from contamination. It has been validated against a bacterial and viral challenge over 24 hours.


This product excels as a bacterial and viral filter and can be used as such in respiratory system providing the maximum period of use in 24 hours, corresponding to the usual correct procedures concerning breathing filters.

**Please note:** This product is contraindicated for use at the patient connection of the breathing system.

### Respiratory system protection



### Oxygen concentrator protection



Code	1790000
Qty.	50
Filtration efficiency	>99.9999%
Resistance to flow at 30L/min	0.9cm H <sub>2</sub> O
Resistance to flow at 60L/min	1.9cm H <sub>2</sub> O
Compressible volume	120ml
Weight	56g
Connectors	22ID – 22OD/15ID
Minimum Tidal Volume	>360ml

[Make an inquiry](#)[Find out more](#)

US Lit No. 82.5 • 03.23