



Intersurgical Venturi valves – Minimum flow rates


Product Code: 0024000
Venturi valve 24% oxygen
Blue
2 L/Min



Product Code: 0028000
Venturi valve 28% oxygen
White
4 L/Min




Product Code: 0031000
Venturi valve 31% oxygen
Orange
6 L/Min




Product Code: 0035000
Venturi valve 35% oxygen
Yellow
8 L/Min




Product Code: 0040000
Venturi valve 40% oxygen
Red
10 L/Min



Product Code: 0070000
Venturi valve 50% oxygen
Pink
15 L/Min



Product Code: 0060000
Venturi valve 60% oxygen
Green
15 L/Min



Quality, innovation and choice



Interact with us
www.intersurgical.co.uk

ventpostint • 12_20