



Intersurgical i-flo™ High flow nasal cannula



Critical Care ▪ i-flo - High flow nasal cannula



Intersurgical i-flo™ - High flow nasal cannula

The i-flo cannula is Intersurgical's single use patient interface designed for the delivery of High Flow Oxygen Therapy (HFOT). The i-flo cannula is available in three sizes and is for use with adult patients in a hospital environment.

HFOT has become a key mode of respiratory support in patients with hypoxemic respiratory failure and is considered to offer a number of benefits versus more conventional oxygen therapies. These benefits have been shown to include:

- Improved oxygenation
- Reduced work of breathing
- Improved secretion clearance
- Improved patient tolerance
- Avoidance of intubation^[1]

The i-flo cannula complements Intersurgical's extensive range of Critical Care breathing systems and interfaces and the unobtrusive design has been developed with a focus on performance and patient comfort.


The i-flo includes a number of unique features and benefits to help ensure the best possible delivery of high flow oxygen therapy.



[1] Masaji Nishimura. High-Flow Nasal Cannula Oxygen Therapy in Adults: Physiological Benefits, Indication, Clinical Benefits, and Adverse Effects. Resp.Care 2016 Apr;61(4):529-41



Adaptive and secure



The i-flo's unique, 180° swivel elbow enables easy and convenient positioning of bedside equipment.



The i-flo clip and lanyard tube management system securely supports the breathing system tubing and allows the patient to move freely without concerns of cannula dislocation or discomfort.

Comfortable and supportive



An anatomical face-mount design directs the head strap above the patient's ears, helping to minimise rubbing and discomfort.

The flexible, semi-rigid face-mount also ensures that strap pressure is distributed evenly over the patient's cheeks. This increases cannula stability and support and helps to prevent pressure hot-spots from developing.





Intersurgical i-flo™ cannula

Code	Description	Box Qty.
7300000	i-flo high flow nasal cannula, small adult	7
7301000	i-flo high flow nasal cannula, medium adult	7
7302000	i-flo high flow nasal cannula, large adult	7



[Make an enquiry](#)

▶ i-flo setup video is available at www.intersurgical.nl




Single limb breathing systems


Code	Description	Box Qty.
2060310	22mm Flextube single heated wire limb with auto-fill chamber and extra limb. ≥ 1.6 m	7



Humidification chamber



Tube management clip



Intersurgical Benelux B.V. Vluchttoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print

Save energy and paper.

If you must print this information sheet please print it double sided.

IS1.53 NL • Issue 1 12.20