








Intersurgical Guedel Airways

New ISO 5364:2016 Specified Colours

Product Code	Intersurgical Size	ISO Size	Previous Intersurgical Colour		New ISO Compliant Colour	
1000035 (S*)	000	3.5	Pink		Light Green	NEW
1100050 (S*)	00	5.0	Blue		Blue	
1110055 (S*)	0	5.5	Grey		Grey	
1111065 (S*)	1	6.5	White		Brown	NEW
1111570 (S*)	1.5	7.0	Yellow		White	NEW
1112080 (S*)	2	8.0	Green		Green	
1113090 (S*)	3	9.0	Orange		Yellow	NEW
1114100 (S*)	4	10.0	Red		Red	
1115120 (S*)	5	12.0	Purple		Purple	

This colour change effects all of Intersurgical's Guedel Airway products, including bulk-packed options.
(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S



ISO Size 3.5



ISO Size 6.5



ISO Size 7.0



ISO Size 9.0

The four sizes of Intersurgical Guedel airways that have had their colours changed to conform to the updated standard.

How to distinguish between old and new colours

If you are unsure whether your Guedel is compliant with the new ISO colour guidance, please compare the colour and ISO size printed on the flange of the Guedel with the table above.

For more information, please visit www.intersurgical.com/info/guedelairway



Quality, innovation and choice



Interact with us

www.intersurgical.com