

Hydro-Guard® Mini

Low Volume Pleated membrane filter - Sterile option



INTERSURGICAL®
COMPLETE RESPIRATORY SYSTEMS

Quality, Innovation and Choice

Hydro-Guard® Mini - sterile option

Intersurgical are pleased to announce that we are now able to provide the Hydro-Guard Mini low volume pleated membrane filter as sterile. This product provides superior bacterial and viral efficiency offering the ultimate in patient protection.

Production of the Hydro-Guard Mini

The Hydro-Guard Mini is manufactured using a two-shot injection moulding process allowing us to create in a single step the integral monitoring port cap and coloured casing.

This eliminates the need for a separate label or port cap and increases the inherent safety of the product in use.

The Hydro-Guard Mini is ideal for use on adult patients in low flow anaesthesia as a heat and moisture exchanging filter and is also suitable for use in intensive care as a filter only product.



Hydro-Guard® Mini - A versatile low volume pleated membrane breathing filter for use in anaesthesia as an HMEF, moisture return 23mg H₂O/L at: VT 500ml, or in ITU as a filter only

Bacterial Filtration efficiency	Viral Filtration efficiency	Resistance to flow filter only		Compressible volume filter only	Weight filter only	Minimum tidal volume filter only
		30L/min	60L/min			
>99.99999%	>99.999%	1.3cm H ₂ O	2.7cm H ₂ O	63ml	30g	200ml

1745000S	125	1745011S	50
<p>22F/15M 22M/15F</p> <p>Luer lock port</p>		<p>22F/15M Superset catheter mount 22M/15F</p> <p>Luer lock port</p>	
Hydro-Guard Mini + luer lock port		Hydro-Guard Mini + luer lock port and Superset catheter mount	

IS12.18 Issue 2 GJ 10.07

United Kingdom (Head Office):
Crane House, Molly Millars Lane,
Wokingham, Berkshire,
RG41 2RZ

T: +44 (0)118 9656 300
F: +44 (0)118 9656 356
info@intersurgical.com
www.intersurgical.com

France
info@intersurgical.fr

Deutschland
info@intersurgical.de

España
info@intersurgical-es.com

Portugal
info@intersurgical.pt

Nederland
info@intersurgical.nl

Česká Republika
info@intersurgical.cz

Philippines
info@intersurgical.ph

Россия
info@intersurgical.ru

Lietuva
info@intersurgical.lt

South Africa
info@intersurgical.co.za

Japan
info@intersurgical.co.za



Certified to ISO 9001 and ISO 14001