

STARMed CaStar Up

Openable hood for CPAP therapy



STARMed CaStar Up – openable CPAP hood

StarMed CaStar Up

Openable CaStar Up is a device recommended for the treatment of respiratory failure:

- Cardiogenic pulmonary oedema
- Non cardiogenic pulmonary oedema (ALI-ARDS)
- Thoracic trauma
- Atelectasis
- Post operative hypoxemia
- Asthma
- COPD exacerbation
- Pre and post extubation respiratory support

Advantages

- Allows rapid access to the patient's face during nursing procedures or in case of emergency
- One operator can easily remove the upper part of the hood and have full access to the patient
- Probes and catheters can remain in place, on the patient, connected to the lower part of the hood
- Well tolerated for extended periods
- Prevents the risk of pressure sores from face masks

Technical Features

- 22M/22F connectors compliant to ISO standard
- Fitted with integrated pressure gauge to check the presence of pressure within the hood
- Internal volume 15 litres – weight 610g
- Comfortable underarm fastening straps
- Other fastening accessories are available if required
- Two new larger access ports on the collar for probes or catheters of 3.5-7.0mm Ø
- Patient access port with screw cap. Fitted with bidirectional anti-suffocation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated



Features and benefits



CE certified product

Transparent

For ease of patient monitoring

Integrated pressure manometer

Anti-suffocation valve

Easy removal for rapid patient access

New large dimension accesses on the collar

Ideal for single and multiple lumen probes or catheters

CaStar Up CPAP hood with anti-suffocation valve and patient access port, ear plugs and measuring tape

Code	Size	Neck Circumference	Box qty.
CA011XXL/2W	Extra Extra Large	≥50cm	5
CA011XL/2W	Extra Large	45-52cm	5
CA011L/2W	Large	40-47cm	5
CA011M/2W	Medium	34-41cm	5
CA011S/2W	Small	27-34cm	5
CA011XS/2W	Extra Small	17-27cm	5

Single packed. Minimum purchasing package: 5 units. Custom-packs are supplied on specific request of the customer.

IS1.31_UK • Issue 1_05.13



Crane House, Molly Millars Lane,
Wokingham, Berkshire, RG41 2RZ
T. +44 (0)118 9656 300
F. +44 (0)118 9656 356

UK
info@intersurgical.co.uk
Italy
info@intersurgical.it

España
info@intersurgical-es.com
Portugal
info@intersurgical.pt

Lietuva
info@intersurgical.lt
Россия
info@intersurgical.ru

France
info@intersurgical.fr
South Africa
info@intersurgical.co.za

Japan
info@intersurgical.co.jp
Nederland
info@intersurgical.nl



Manufactured by: Intersurgical S.p.A.