



# Silente™ adjustable venturi valve

Low operating noise



Oxygen and Aerosol Therapy ▪ Fixed Oxygen Concentration (High Flow)



## Silente™ adjustable venturi valve

The Silente has been designed to operate at a very low noise level. The large diameter body provides space for clear, widely spaced calibration and incorporates a 'pull-twist' action to prevent accidental adjustment.

The twin-jet design offers the user seven oxygen concentrations. The increased body diameter provides space for clear, widely spaced calibration and incorporates a 'pull-twist' action preventing accidental adjustment. The one-piece design eliminates the need to use more than one valve for all concentrations and there are no separate pieces to get lost or mislaid.

The concentration delivered is unaffected by the patient's breathing pattern provided that the optimum flow indicated on the valve is set on the oxygen flow meter.

### Ease of use

The Intersurgical adjustable venturi valve is simple to use and has the following features and benefits:

- The appropriate flow range is selected to suit the desired concentration
- The outer ring with its 'pull-twist' action selects the prescribed oxygen concentration
- A flow in L/min is printed above each concentration indicating the optimum flow meter setting

### Precision

The Silente delivers precise oxygen concentrations. This accuracy is independent of the oxygen flow meter setting and is achieved using just one valve.

### Incurved nose seal

Conforms to different nose shapes designed to prevent oxygen entering patient's eyes

### Large, well spaced calibration settings

Easy to read, minimizes errors

### 'Pull-Twist' controls

Prevent inadvertent adjustment

### Extremely quiet

Improved patient environment

### Soft face seals

Increases patient comfort

### Elastic can be positioned under or over the ears

Below ear position eliminates trauma to top of ears

### One piece design

No parts to change or lose

**A choice of "under chin" or "on chin" positions**  
Provides a better fit on a wider range of patient face shape



Make an inquiry

Code	Description	Tube length	Box Qty.
1137015	Silente, Intersurgical EcoLite™, adult, adjustable venturi mask kit with tube	1.8m	30
1138015	Silente, Intersurgical EcoLite, pediatric, adjustable venturi mask kit with tube	1.8m	40
1188015	Intersurgical EcoLite, adult aerosol mask		50
1198015	Intersurgical EcoLite, pediatric, aerosol mask		42
1139030	Silente adjustable venturi, valve only		50

Lower environmental impact product

US Lit No. 86.5 • Issue date 05.19