



Right-angled water traps

Self-sealing water traps



Critical Care ▪ Accessories

Right-angled water traps

Self-sealing water traps

Intersurgical's right-angled water traps for the collection of condensation compliments our existing range of products, the angled water trap is available in two sizes.

Medium volume and large volume, both are ideal for use in applications where a 90° attachments may be required.

The self-sealing facility allow the water traps to be emptied without interrupting gas flow.

Both water traps feature a 22ID machine side connection to a 22OD breathing system side connection, and can be attached directly to a ventilator or to an anesthetic machine.



Self-sealing facility

Allow the water traps to be emptied without interrupting flow



Code	Description	Qty.
1916000	Water trap, self-sealing, right-angled, medium volume	55
1915000	Water trap, self-sealing, right-angled, large volume	40

Make an inquiry