



VariFit™ – NIV mask



Critical Care ▪ Patient Interfaces



VariFit™ featuring AIR^ogel® technology

The VariFit is a single patient use hospital NIV mask featuring a cushion with AIR^ogel® technology for enhanced patient comfort. Colour coded elbows help the clinician to choose the correct mask for bilevel/CPAP using a single limb breathing system, or non-invasive ventilation using a dual limb breathing system. The breathable, quick release headgear is designed to fit easily to the patient and improve patient comfort.



AIR^ogel® Maximum comfort, minimum weight.
AIR^ogel® minimises pressure points and skin irritation

Advanced Cushion Technology™

The ergonomic cushion technology forms an effective seal for a comfortable fit

Active Headgear Connector™

Moves with the patient for a reliable seal

360° rotating elbow

For optimal tube positioning and easier mobility

Flexible Spacebar™

The unique bendable design provides superior stability, maintaining a seal without contacting the skin

Vented Mask (only 2260000 to 2262000)

Vented mask option with an anti-asphyxiation valve for use with CPAP or bilevel devices

Oxygen Port

Allows easy delivery of supplementary oxygen

Colour Coded Elbows

To help the clinician choose the correct option

Anti-asphyxiation Valve

Prevents rebreathing in the case of inadequate ventilation

Vented Mask with an anti-asphyxiation valve for use with a single limb breathing system that does not have a CO₂ leak port

Clear elbow ○



Non-vented Mask with an anti-asphyxiation valve for use with single limb breathing systems that have a CO₂ leak port

Clear and blue elbow ◐



Non-vented Mask for use with dual limb breathing systems

Blue elbow ●



Code	Description	Colour coded elbow	Box Qty.
For use with a single limb breathing system that does not have a CO₂ leak port			
2260000	VariFit NIV mask, vented with anti-asphyxiation valve, small adult	○	8
2261000	VariFit NIV mask, vented with anti-asphyxiation valve, medium adult	○	8
2262000	VariFit NIV mask, vented with anti-asphyxiation valve, large adult	○	8
For use with a single limb breathing system that has a CO₂ leak port			
2263000	VariFit NIV mask, non-vented with anti-asphyxiation valve, small adult	◐	8
2264000	VariFit NIV mask, non-vented with anti-asphyxiation valve, medium adult	◐	8
2269000	VariFit NIV mask, non-vented with anti-asphyxiation valve, large adult	◐	8
For use with dual limb breathing systems			
2265000	VariFit NIV mask, non-vented, small adult	●	8
2266000	VariFit NIV mask, non-vented, medium adult	●	8
2267000	VariFit NIV mask, non-vented, large adult	●	8

Make an enquiry

AIRgel is a registered trademark of Sleepnet Corporation. Flexible Spacebar, Advanced Cushion Technology and Active Headgear Connector are trademarks of Sleepnet Corporation.

IS1.41_UK • Issue 3 06.16



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com

UK • Ireland • France • Germany • Spain • Portugal • Italy
Netherlands • Sweden • Denmark • Lithuania • Russia
Czech Republic • Turkey • South Africa • China • Japan
Taiwan • Philippines • USA • Colombia • Australia



Please think before you print
Save energy and paper.
If you must print this e-mail
please print it double sided

Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485